


Canon



Aplio a550

Advanced. Integrated.
Seamless.



Aplio a550



Advanced. Integrated. Seamless.

Aplio a550 integrates industry-leading imaging technologies, advanced applications and intuitive controls for the busy clinician. Designed to enhance productivity and throughput while maximizing clinical confidence, Aplio a550 is the ideal system for routine imaging departments. Its extensive selection of advanced applications provides confidence your system will meet a wide range of clinical needs.

Boost your clinical confidence

Aplio's powerful imaging technologies provide you with better image quality with reduced clutter, strengthened signal and improved visualization. Aplio a-series' unique aBeam architecture provides the capability to ensure that all of Aplio's unique imaging technologies work together seamlessly for greater uniformity across all applications.



Precision+ offers outstandingly smooth images with sharpened outline of lesions, enhanced image uniformity and reduced clutter.




ApliPure+ compounding delivers increased imaging contrast and reduced speckle noise to improve visualization.




Better diagnostics starts here

Designed to increase efficiency, the system's lightweight transducers feature outstanding clinical versatility, ergonomic shapes and thin, super-flexible cables. Aplio a550 is compatible with a wide range of transducers from across the Aplio product range, ensuring high productivity while helping reduce cost for specialty probes.




Differential Tissue Harmonics provides harmonic images of unsurpassed spatial resolution, alongside greatly enhanced penetration.

Aplio's wideband transducer and signal processing technology delivers outstanding sensitivity, penetration and spatial resolution for all Doppler modes.



SMI's level of vascular visualization, combined with high frame rates, advances diagnostic confidence when evaluating lesions, cysts and tumors.





Color-coded SMI allows you to depict flow and greyscale information with high temporal and spatial information simultaneously.



Seeing the unseen


Experience color flow imaging with unmatched detail and definition on Aplio a550. Superb Micro-vascular Imaging (SMI) expands the range of visible blood flow to a level of detail never before seen with diagnostic ultrasound. Exclusive iSMI now also allows you to image larger regions of interest without a reduction in frame rate.




Traditional color Doppler imaging (left) removes clutter from the images by suppressing low-velocity components, resulting in a loss of flow in tiny vessels. SMI (right) separates flow from overlaying tissue motion effectively, while preserving even the subtlest low-flow components with unmatched detail and definition.

Increase your confidence, expand your capability

Aplio a-series' suite of advanced imaging and quantification functions provides the metrics to help you quickly diagnose with confidence. Combined with tools for early detection and reliable characterization of lesions can help you optimize your patients' clinical pathway.



Canon Medical Systems' shear wave technology provides a quantitative measure and realtime display of tissue elasticity in a variety of clinical settings ranging from abdominal to small parts examinations.




Smart Maps help you visualize and quantify shear wave propagation in realtime. Aplio's unique propagation map is a powerful and intuitive tool to visually assess the quality of an elastogram.

The system's comprehensive strain elastography suite with raw data functionality assists you in localizing and assessing palpable masses with high accuracy, sensitivity and reproducibility.


Aplio's comprehensive CEUS imaging and quantification package allows you to assess perfusion dynamics in a wide range of clinical settings, including an ample variety of specialized exams.

Navigate with ease, treat with confidence

Aplo provides a wide range of tools for advanced imaging and interventions. Dedicated transducers and an abundance of imaging and navigation tools help you enhance confidence and accuracy during interventional procedures and their follow-up.




Smart Fusion merges realtime ultrasound with previously acquired CT, MR or ultrasound data, allowing you to identify and compare lesions with ease and to navigate complex anatomy confidently.



For comprehensive evaluation, Smart Fusion allows you to work in multiple imaging modes, including color Doppler and CEUS. The concise quad display shows the live ultrasound image in sync with multiple views of the pre-loaded data.

Aplio a-series is compatible with a variety of needle guides with multi-angle or free angulation capability, either using brackets or directly mounted on the transducer to ensure easy handling with high precision and minimal blind zone.




Merging modalities to improve visualization

On Aplio, merging two modalities for synchronized display is a simple and quick two-step process. Intelligent target and marker points facilitate navigation in the region of interest.



Better intercostal access

Aplio's thin convex transducers are ideally suited for intercostal scanning. The new biopsy attachment with minimized blind area and selectable puncture angle facilitates optimal puncture conditions for each patient.



Exceptional detail for a more precise diagnosis


Both the busy clinician and the patient can benefit from high-resolution 2D imaging and volumetric ultrasound. Aplio's comprehensive volume imaging suite extends your diagnostic capabilities into the next dimension of imaging with extraordinary image quality and uncompromised workflow.



Doppler Luminance provides a homogeneous, easy-to-interpret color display with high accuracy and rich detail, even in the smallest vessels. Doppler Luminance offers high frame rates while maintaining full B-mode image quality.



Luminance offers natural-looking 3D renderings of high quality and definition, providing strong visual feedback on depth and detail from the first trimester onwards.





A wide range of MultiView options provides high-resolution cross sections, helping you to better understand anatomical relationships or the extent of a given structure.

Rt Follicle(3D)	
1	2.61 ml
2	1.08 ml
3	0.71 ml
4	0.61 ml
5	0.15 ml
6	0.15 ml
7	0.13 ml
8	0.08 ml
9	0.04 ml
10	0.03 ml


Accurate quantification, regional myocardial function

Functional assessment is at the heart of cardiovascular imaging. By providing valuable additional information in easy-to-understand visual, parametric or quantitative formats, Aplio's advanced clinical functions can help you get your diagnostic answer faster and more reliably.




Aplio's clinical protocol assistant provides a reliable method to ensure that complex exams are done consistently patient after patient. Once activated, a clear, easy-to-read menu will guide you step by step through the entire exam.

Aplio provides you with high frame rate Tissue Doppler images and Pulsed-Wave-TDI traces for a precise timing of cardiac events in both visual and quantitative formats.



Supporting standard and user-defined protocols for both physical and pharmacological stress, Aplo offers a comprehensive package for fast and accurate wall motion assessment.



Aplo's advanced Wall Motion Tracking technology provides immediate visual and quantitative access to global and regional myocardial wall motion dynamics.



Aplio makes your work flow

Aplio provides a host of intelligent workflow support and automation tools, helping you to achieve rapid results with consistent high quality regardless of the patient condition.




Designed with our users in mind

Smaller and lighter, Aplio a550 is easy to maneuver. With over 36 cm panel height adjustment, lateral slide and a fully articulating monitor arm, Aplio a550 helps you to optimally adjust the console to virtually any scanning position.

Switch to auto-pilot

Aplo's context-sensitive user interface is designed to make your imaging task simpler and quicker. While automated settings can deal with routine clinical needs, you always retain control over all imaging parameters when needed.



The mode-sensitive on-screen navigation for the central trackball boosts your workflow and efficiency. By visually guiding you through the exam, it allows you to adapt and operate the system within a few minutes.

Access all areas


Aplio's large, tablet-style touch screen with three interactive zones allows you to quickly browse and select the desired function, while the rest of the display remains unchanged.



Realtime QuickScan allows you to achieve greater consistency in your exams by ensuring that superb image quality is the benchmark.



Aplio's embedded raw data architecture allows you to optimize, review, analyze and report clinical data anytime with no loss of functionality.



A range of automated measurement and analysis tools helps you increase accuracy, consistency and speed of your exams.

Aplio *a-series*



CANON MEDICAL SYSTEMS EUROPE B.V.

<https://eu.medical.canon>

©Canon Medical Systems Corporation 2019. All rights reserved.
Design and specifications subject to change without notice.
Model number: CUS-AA550
MCAUS0303ECA 2019-04 CMSE/Printed in Europe

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485.
Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Aplio, ApliPure, Dynamic Flow and Made for Life are trademarks of Canon Medical Systems Corporation.

Disclaimer: Some features presented in this brochure may not be commercially available on all systems shown or may require the purchase of additional options. Please contact your local Canon Medical Systems representative for details.

Made For life